

Mouse Pan B Cell Isolation Kit (Column Free,Negative) Instructions

Product Information

Product Name	Model	Product Components	Specification
Mouse Pan B Cell Isolation Kit-100T (Column Free,Negative)	F1501-10	Mouse Pan B Cell Isolation Cocktail	1 mL
		Streptavidin MicroBeads	1 mL
Mouse Pan B Cell Isolation Kit-20T (Column Free,Negative)	F1501-10T	Mouse Pan B Cell Isolation Cocktail	200 µL
		Streptavidin MicroBeads	200 µL

Product Description

The Mouse Pan B Cell Isolation Kit (Column Free,Negative) (the “Kit”) is designed for labeling the non-Pan B cells using isolation cocktail and streptavidin microbeads, and removing the labeled cells via the RWD magnetic pole, achieving column-free cell isolation.

The Kit enables rapid, easy, and high-yield isolation of Pan B cells from spleen tissue of mice. The isolated cells are free of antibody and microbeads labeling, can be directly used for downstream applications including flow cytometry and cell culture.

Test Capacity

- ✧ F1501-10 can process up to 1*10⁹ of total cells.
- ✧ F1501-10T can process up to 2*10⁸ of total cells.

Storage & Transportation

- ✧ The Kit should be transported at 2 ~ 8°C.
- ✧ The Kit should be stored at 2 ~ 8°C. Protect from light. Do not freeze.
- ✧ The Kit is valid for 12 months from the date of manufacture.

Consumables & Instruments

Reagents	Consumables	Instruments
Buffer	Magnetic Pole (RWD: #CFMS-S1)	Refrigerated Microcentrifuge (RWD)
High Efficiency Debris Removal Kit (RWD)	5 mL Flow Cytometry Tube	Fluorescent Automated Cell Counter (RWD)
	40 µm Cell Strainer	Single Cell Suspension Dissociator (RWD)

Operation

Preparation before Experiment

- Preparation of Buffer: Filter PBS (pH 7.2 ~ 7.4, containing 0.5% BSA and 2 mM EDTA) through a 0.22 µm membrane under 2 ~ 8°C for sterilization and store at 2 ~ 8°C.
⚠ Note: PBS containing Ca²⁺or Mg²⁺ or biotin is not recommended.

Sample Preparation

- For peripheral blood samples, isolate peripheral blood mononuclear cells (PBMC) by density gradient centrifugation;

For solid tissue samples (e.g., spleen and lymph node), before cell isolation (enrichment), dissociate the tissue into a single cell suspension using the single cell suspension dissociator.

- Filter the cell suspension through the 40 µm cell strainer to remove debris and cell clumps.
⚠ Note: It is not recommended to perform red blood cell lysis, which may compromise the cell viability and surface antigen integrity. Since the dead cells may have non-specific binding with antibody or microbeads, it is recommended to perform density gradient centrifugation or use High Efficiency Debris Removal Kit (RWD: #DHDR-5006) to remove dead cells.

Cell Labeling

- ⚠ Note: The recommended reagent volumes below can process 1*10⁷ of total cells. If the cell count is less than 1*10⁷, use the same volume as for 1*10⁷ cells. If the cell count exceeds 1*10⁷, increase the reagent volumes proportionally.
- Cell Concentration Adjustment:** Adjust the cell concentration to 1*10⁸ cells/mL.
 - Antibody Labeling:** Add 100 µL of sample and 10 µL of mouse Pan B cell isolation cocktail to the 5 mL flow cytometry tube. Mix them evenly and incubate at room temperature for 5 min.
⚠ Note: Add no more than 1 mL of sample per tube (cell count not exceeding 1*10⁸).
 - MicroBeads Labeling:** Add 10 µL of streptavidin microbeads to the sample.
⚠ Note: For optimal labeling efficiency, gently mix the suspension 2 ~ 3 times during incubation.
 - Volume Standardization:** Add the buffer to adjust the total sample volume to 2.5 mL. Mix evenly and incubate at room temperature for 2 min.

Cell Isolation

- Insert the flow cytometry tube containing sample to the magnetic pole for magnetic capture at room temperature for 3 min.
- Take out the magnetic pole and pour the sample to a new 5 mL flow cytometry tube. After pouring, pull out the tube from the magnetic pole.
⚠ Note: During pouring, invert the magnetic pole with the flow cytometry tube for 2 ~ 3 s. Do not pull out the tube from the pole until the pouring is complete. After pouring, directly discard the remaining liquid in the tube.

Results

Isolate the Pan B cells from the spleen cells of C57 mice using the Mouse Pan B Cell Isolation Kit (Column Free,Negative) in combination with the RWD magnetic pole (RWD: #CFMS-S1). The isolated cells are labeled with PE anti-mouse CD19 Antibody (Clone: #6D5) and APC anti-mouse CD138 (Syndecan-1) Antibody (Clone: #281-2) followed by flow cytometric analysis. The percentages of Pan B cells (CD19+CD138-) before and after isolation are 63.81% and 98.44% respectively.




Figure (a) Before vs. After Isolation

Precautions

- (1) The Kit is for scientific research only, not intended for any diagnostic or therapeutic use.
- (2) The Kit is valid for 12 months and RWD shall not guarantee the validity of expired products.
- (3) Ensure that all steps are performed under sterile conditions.
- (4) Do not mix kits of different models/lots.
- (5) To improve the purity of magnetic-labeled cells, add 2.4 mL of buffer to the new flow cytometry tube in the last step of **Cell Isolation**. Mix evenly and insert the tube containing sample to the magnetic pole for magnetic capture at room temperature for 2 min. Then, repeat step (2) of **Cell Isolation**.

© 2026 RWD Life Science Co., Ltd. All Rights Reserved.

RWD Life Science Co., Ltd.

Add: 10410 Corporate Drive, Sugar Land, TX 77478, USA

Add: (Floor 9, 19&20 Building 7A, Floor 9 Building 7D) Room 1901, Building 7A, International Innovation Valley, Dashi 1st Road, Xili Community, Nanshan District, Shenzhen 518000, Guangdong, P. R. China

Web: www.rwdstco.com

E-Mail: service@rwdls.com

Tel: 0086-755-86111281

001-858-900-6602(USA)