

RWD VETERINARY PRODUCT CATALOG

RWD Life Science Co.,Ltd

Add: 19-20/F, Building 7A, Shenzhen International Innovation Valley, Dashi 1stRoad, Nanshan District, Shenzhen, Guangdong, P.R. China E-mail: rwd@rwdstco.com

RWD Life Science Inc.

Tel: (858)900-5879 Support: service@rwdls.com Web: www.rwdstco.com

Copyright ©RWD Life Science Co., LTD. 2020. All rights reserved.

No part of this document maybe reproduced or transmitted in any form or by any means without prior written consent of RWD Life Science Co., LTD.

Trademark Notice

RWD, RWD | ******, and RWD are trademarks or registered trademarks of RWD Life Science Co., LTD.

Other trademarks, product, service and company names mentioned are the property of their respective owners.

General Disclaimer

The information in this document may contain predictive statements including, without limitation, statements regarding the future financial and operating results, future product portfolio, new technology, etc. There are a number of factors that could cause actual results and developments to differ materially from those expressed or implied in the predictive statements. Therefore, such information is provided for reference purpose only and constitutes neither an offer nor an acceptance.

RWD may change the information at any time without notice.

www.rwdstco.com



Introduction

Founded in 2002, RWD Life Science is specialized in Animal Healthcare, Life Science and Medical Devices. Headquartered in Shenzhen, our products and services have been exported to over 100 countries and regions.

While dedicated to improving the quality of life, we are also firmly committed to being a leading global science enterprise. Oriented by customer demands and motivated by excellent talents, RWD is always working together with labs & universities, and investing into R&D to achieve the innovation in our products, technologies and service, and to provide our customers with higher- quality solutions.

Being based in China and serving the world, RWD has not only achieved a strong reputation with advantages in the field of veterinary inhalation anesthesia system, lab research and surgery instrument, but also gained long-term trust and support from customers all over the world in the past 18 years.

Looking into the future, we shall adhere to our core values "Quality, Integrity, Responsibility", and also contribute to the development of pre-clinic and medical clinic applications in the international market.





Content

| 01/ | D-Pro Veterinary Dental Unit Veterinary Dental X-ray | 01 |
|-----|--|----|
| 02/ | Anesthesia Machine R620-S1 Veterinary Anesthesia Machine R640-S1 Economical Veterinary Anesthesia Machine RSA Small Pet Anesthesia Machine R640 Flexible Veterinary Anesthesia Machine-Wall Mount R620-MRI MR Conditional Animal Anesthesia Machine R630 Dual-drug Veterinary Anesthesia Machine R610 Veterinary Anesthesia Machine-Tabletop | 03 |
| 03/ | Ventilator R419 Intelligent Animal Ventilator | 10 |
| 04/ | Monitor RM300 / RM400S / RM400M RM-C1 Veterinary Capnograph R-CU500 The Smart Animal ICU | 11 |
| 05/ | Gas Supply ROC-5A Veterinary Oxygen Concentrator Anesthesia Air Pump | 13 |
| 06/ | Waste Gas Scavenging Automatic Weighing Meter Gas Filter Canister | 14 |
| 07/ | Accessories Vaporizer Anesthesia Induction Chambers Shelf Laryngoscopes Toggle Switch Regulator | 15 |
| 08/ | Disposable Mask Breathing Bags Breathing Circuits Endotracheal Tubes Air Hose Urinary Catheter | 17 |
| 09/ | Animal Surgical Instruments | 21 |
| 10/ | One-Stop Solution | 23 |



D-Pro Veterinary Dental Unit

D-pro veterinary dental unit provides a variety of functions for veterinarians, such as grinding, polishing, cutting and cleaning, and makes daily animal oral health work easy and efficient.





Veterinary Ultrasonic Scaler (optional)

- Water and electricity separation design to ensure the safety of the equipment.
- 360 ° hand piece, light in weight, reduces hand fatigue effectively.
- High-speed LED hand piece provides bright operative vision.
- Low-speed hand piece helps veterinarians to further clean the plaque and inhibit the growth of the plaque.
- High quality oil-free air compressor to avoid the trouble of daily maintenance
- 20L large capacity air tank ensures pressure balance and long-term supply, and reduces the frequency of air charging.
- Magnetic holder to prevent fine dental tools accidental loss
- The whole machine weighs only 46kg, flexible movement with super smooth casters

Order Information

| Cat No. | Description |
|---------|------------------------|
| D-pro | Veterinary Dental Unit |

Included

| Description | Amount | Function | Remarks |
|------------------------------|---|---|----------|
| Host | 1pcs | | |
| High Speed Handpiece | 1pcs | Mainly for Cutting | |
| Low Speed Handpiece | 1pcs | Polishing | |
| 3 Way Air/Water Syringe | 1pcs | Oral Disinfection and Washing | |
| Suction Handpiece | 1pcs | Take away the liquid in the mouth | Included |
| Disposable Suction Tip | 1Set (25pcs) | Used with Suction Device | |
| Burs | 1Set (Tapered fissure burs*5+Spherical burs*5) | Used with High Speed Handpiece | |
| Disposable Polishing Head | 1Set (10pcs) | Used with Low Speed Handpiece | |
| Polishing Paste | 1Set (10pcs) | Used with Polishing Head | |
| Lubricant for Handpiece | 1bottle | 30ml ,For daily maintenance of high / low speed Handpiece | Optional |
| Veterinary Ultrasonic Scaler | 1Set (6 Scaling Tips Included) | For Animal Teeth Scaling | |

Veterinary Dental X-ray

Used for the diagnosis of pet teeth. Veterinarians make the treatment plan according to the X-ray image.



- Using the international advanced technology, efficient integrated design
- Microcomputer intelligent control, not only remote controlled exposure, but more powerful function of low voltage alarm and high voltage protection
- Micro focus technology, much more clear image
- The dual focus configuration makes the image clearer and the diagnosis more accurate
- Can use lightroom dental film, imaging in one minute
- Intraoral digital imaging system is optional

Order Information

| Model | Description |
|-------|--------------------------------|
| D-X8V | Mobile Veterinary Dental X-ray |

Included

| Description | Amount | Remarks |
|-------------------------------------|--------|----------|
| Mobile Veterinary Dental X-ray Host | 1 | |
| Exposure Switch | 1 | |
| Remote Control | 1 | |
| Dental digital imaging system | 1 | Optional |



▲ D-X2V

- Small size, light weight, portable
- With memory storage function, only one exposure can permanently memorize the exposure parameters
- 0.4mm * 0.4mm dual focus configuration, which makes the image clearer and the diagnosis more accurate
- High-frequency power supply Low radiation Good-looking appearance
- Can use lightroom dental film, imaging in one minute

Order Information

| Model | Description |
|-------|----------------------------------|
| D-X2V | Portable Veterinary Dental X-ray |

Included

| Description | Amount | Remarks |
|---------------------------------------|--------|-------------|
| Portable Veterinary Dental X-ray Host | 1 | |
| Mobile Scaffold | 1 | Optional |
| Dental Digital Imaging System | 1 | - Coptional |

01 RWD Life Science | www.rwdstco.com www.rwdstco.com | RWD Life Science | 02

AMD



R620-S1 Veterinary Anesthesia Machine

R620-S1 is applicable to clinical anesthesia of dog, cats, pigs, monkeys, rodent, reptiles, birds etc. With minimal operation design, vet shall feel easy, safe and reliable to perform surgery on animals.









Circuit switch

Rail-mounted CO₂ absorber

APL valve

- High accuracy vaporizer
- APL valve for quick-occlusion
- Rail-mounted CO₂ absorber, so easy to install and disassemble
- One switch to shift between rebreathing and non-rebreathing circuits
- Shelf for placing monitor and ventilator
- Good to hang E-type oxygen cylinder, and oxygen concentrator
- Good waste gas scavenging system
- Complete function with both breathing and non-rebreathing circuits
- Five-foot base/ Four-foot base/Wall-mount Holder is available

Order Information

| Cat No. | Description |
|---------|-------------------------------|
| R620-S1 | Veterinary Anesthesia Machine |

Included

| Cat No. | Description | Remarks |
|--------------|---|------------------------------------|
| R620-S1-A | R620-S1 Series Mounting Assembly with ${\rm O_2}$ Flowmeter (0.1-4LPM)_EN | |
| R580S | Vaporizer for Isoflurane- Pour Fill- Cagemount | Can be changed to other vaporizer. |
| R580-IO | Vaporizer Endcaps (1 pair) | |
| R510-31 | Gas Filter Canister, Large | |
| R600-ET | Cuffed Endotracheal Tubes (1 set) | |
| R600-NL-0.5L | Reusable Non Latex Breathing Bag, 0.5L | |
| R600-NL-1L | Reusable Non Latex Breathing Bag, 1L | |
| R600-NL-2L | Reusable Non Latex Breathing Bag, 2L | |
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |

R640-S1 Flexible Veterinary Anesthesia Machine

R640-S1 is an economical machine designed with good performance but competitive price. Applicable to clinical anesthesia of dogs, cats, pigs, monkeys, rodents, reptiles, birds etc.



- Compact structure & simple operation
- Two colors: red or silver
- Support re-breathing circuit and non-rebreathing circuit (Basic configuration re-breathing circuit 1.5m)
- Pressure gauge Constantly monitoring the circuit pressure
- Oxygen flush button Preventing emergencies
- Red rotary float More eye-catching and readable
- Freewheel with lock More stable
- Five-foot base/ Four-foot base/Wall-mount Holder is available



Order Information

| Cat No. | Description |
|---------|--|
| R640-S1 | Flexible Veterinary Anesthesia Machine |

Included

| Cat No. | Description | Remarks | |
|--------------|---|---------------------------|----------|
| R640-S1-A-S | R640-S1 Series Mounting Assembly (Silver) with O_2 Flowmeter (0.1-4LPM) | Can change to the red one | |
| R580S | Vaporizer for Isoflurane- Pour Fill- Cagemount | Can be changed | |
| R580-IO | Vaporizer Endcaps (1 pair) | | |
| R600-NL-1L | Resuable Non Latex Breathing Bag, 1L | | |
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m | | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | | |
| R640-S1-MS-S | Monitor Shelf (for R640-S1) | Can change to the red one | Ontinual |
| R600-ET | Cuffed Endotracheal Tubes (1 set) | | Optional |
| R510-31 | Gas Filter Canister, Large | | |

RSA Small Pet Anesthesia Machine

RSA inhalation anesthesia machine is specially designed for small pets such as reptiles, birds, cats, and small mammals within 12 kg. It features with dual-flow meter (air & oxygen) to meet the demands of different species.



- 0-4L/min dual-flow meter (air and oxygen)
- Non-rebreathing system specially for animals within 12 kg
- Embedded oxygen flush button to prevent accidentally operating
- High accuracy from vaporizer (tested one by one before delivery)
- Compact structure and easy operation
- Red rotary float---more eye-catching and accurate













Order Information

| Cat No. | Description |
|---------|------------------------------|
| RSA | Small Pet Anesthesia Machine |

Included

| Cat No. | Description | Remarks |
|--------------|---|----------------|
| RSA-A(H)-S | Small Pet Anesthesia Machine Mounting Assembly (Silver) | |
| R580-IO | Vaporizer Endcaps (1 pair) | |
| R583S | Vaporizer for Isoflurane- Easy Fill- Cagemount | Can be changed |
| R-NBC-C(POV) | Modified Mapleson C System (Pop-Off Occlusion Valve) | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |
| R510-31 | Gas Filter Canister, Large | |
| V100-V | Anesthesia Induction Chamber-Small(15*10*10cm) | |
| V101-V | Anesthesia Induction Chamber-Medium (24*12*18cm) | |
| V102-V | Anesthesia Induction Chamber-Large (40*18.5*25cm) | Ontinal |
| R-HS | Oxygen Concentrator Shelf (for H-base) | Optinal |
| R640-S1-MS-S | Monitor Shelf (for R640-S1) | |
| R510-29 | Anesthesia Air Pump (HF), 220V | |
| R510-30 | Anesthesia Air Pump (HF), 110V | |

R640 Flexible Veterinary Anesthesia Machine-Wall Mount

R640 is an excellent choice for Non-operating room anesthesia, like Treatment Room, Imaging Room, Anesthesia Preparation Room and Feline Specialist Clinics. It can effectively improve the treatment rate and the medical experience of the owner and animals.



- Wall mount, ideal if you have limited space.
- Non-rebreathing circuit, ensure fast induction, quicker wake up and recovery.
- 10L flow meter, good for animals less than 50kg.
- Oxygen flush with safety restrictor.

Order Information

| Cat No. | Description |
|---------|---|
| R640 | Flexible Veterinary Anesthesia Machine-Wall Mount |

Included

| Cat No. | Description | Remarks |
|-------------|---|----------------|
| R640-H | R640 Series Mounting Assembly | |
| R582S | Vaporizer for Sevoflurane- Easy Fill- Cagemount | Can be changed |
| R580-IO | Vaporizer Endcaps (1 pair) | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | |
| R600-NL-3L | Reusable Non Latex Breathing Bag, 3L | |
| R600-NL-1L | Reusable Non Latex Breathing Bag, 1L | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |

5 RWD Life Science | www.rwdstco.com www.rwdstco.com | RWD Life Science 06

R620-MRI MR Conditional Animal Anesthesia Machine

R620-MRI has been tested by Shenzhen Institutes of Advanced Technology. And R620-MRI is proved not to produce displacement and not affect the operation and imaging in MRI systems of 3 Tesla and below 3 Tesla.



- Can be used in MRI systems of 3 Tesla and below 3 Tesla.
- · High accuracy vaporizer.
- Good waste gas scavenging system.
- Complete function with both breathing and non-rebreathing circuits.

Order Information

| Cat No. | Description |
|----------|--|
| R620-MRI | MR Conditional Animal Anesthesia Machine |

Included

| Cat No. | Description | Remarks |
|--------------|--|---------|
| R580S-MRI | Vaporizer for Isoflurane- Pour Fill- Cagemount-MRI | |
| R580-IO-MRI | Vaporizer Endcaps (1 pair)- MRI | |
| R510-31 | Gas Filter Canister, Large | |
| R600-ET | Cuffed Endotracheal Tubes (1 set) | |
| R600-NL-0.5L | Reusable Non Latex Breathing Bag, 0.5L | |
| R600-NL-1L | Reusable Non Latex Breathing Bag, 1L | |
| R600-NL-2L | Reusable Non Latex Breathing Bag, 2L | |
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | |
| R-AH-MRI | Air Hose,ID8.0-9/16(SGL),5m,MRI-Compatible | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |
| | | |

R630 Dual-drug Veterinary Anesthesia Machine

R630 was developed according to the requirement of veterinary clinics. Suitable for clinical anesthesia of dogs, cats, pigs, monkeys, rodents, reptiles, birds, etc.







Non-return

H-base(optional)

- Dual vaporizers for isoflurane and sevoflurane
- "Non-return" parallel structure to avoid drug mixing
- Easy Fill to avoid leakage and wrong anesthetic filling
- Excellent user experience with APL vale, circuit switch and rail-mounted CO₂ absorber
- Good waste gas scavenging system
- Complete function with both breathing and non-rebreathing circuits

Order Information

| Cat No. | Description |
|---------|---|
| R630 | Dual-drug Veterinary Anesthesia Machine |

Included

| Cat No. | Description | Remarks |
|--------------|---|------------------|
| R630-A | R630 Series Mounting Assembly with O2 Flowmeter (0.1-4LPM)_EN | |
| R582S | Vaporizer for Sevoflurane- Easy Fill- Cagemount | |
| R583S | Vaporizer for Isoflurane- Easy Fill- Cagemount | |
| R580-IO | Vaporizer Endcaps (1 pair) | 2 pairs included |
| R510-31 | Gas Filter Canister, Large | |
| R600-ET | Cuffed Endotracheal Tubes (1 set) | |
| R600-NL-0.5L | Resuable Non Latex Breathing Bag, 0.5L | |
| R600-NL-1L | Resuable Non Latex Breathing Bag, 1L | |
| R600-NL-2L | Resuable Non Latex Breathing Bag, 2L | |
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |

Ventilator

R419 Intelligent Animal Ventilator

R419 is an intelligent ventilator super easy for vets to operate: pre-setted parameters will automatically pop up by only entering weight of the animal. Support VCV, PIPCV and APNEA control to meet multiple clinical applications.













Adjustable pressure range: 3~20cmH2O

Material: polycarbonate, polypropylene, silica rubber, stainless steel



- Automatic pre-setted parameters by entering animal weight
- 20 sets of parameters can be stored and retrieved by vets
- APNEA mode available, do not interfere patient spontaneous breathing
- Identify and alarm suffocation, monitor the respiratory during wholeprocess
- Touch screen and 75° visual surface design for easy operation













- Tabletop structure and space saving
- High accuracy vaporizer
- High-efficiency gas filter canister
- Support breathing and non-rebreathing circuits
- Oxygen flush for first aid
- Can be extended to mobile anesthesia machine

Order Information

| Cat No. | Description |
|---------|--|
| R610 | Veterinary Anesthesia Machine-Tabletop |

R610 Veterinary Anesthesia Machine- Tabletop

R610 is applicable to clinical anesthesia of dog, cats, pigs, monkeys, rodent, reptiles, birds etc. With minimal operation

design, vet shall feel easy, safe and reliable to perform surgery on animals.

Included

| Cat No. | Description | Remarks |
|--------------|---|----------------|
| R610-A | R610 Series Mounting Assembly (Silver) with O ₂ Flowmeter (0.1-4LPM) | |
| R580S | Vaporizer for Isoflurane- Pour Fill- Cagemount | Can be changed |
| R580-IO | Vaporizer Endcaps (1 pair) | |
| R510-31 | Gas Filter Canister, Large | |
| R600-ET | Cuffed Endotracheal Tubes (1 set) | |
| R600-NL-0.5L | Resuable Non Latex Breathing Bag, 0.5L | |
| R600-NL-1L | Resuable Non Latex Breathing Bag, 1L | |
| R600-NL-2L | Resuable Non Latex Breathing Bag, 2L | |
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m | |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) | |
| R-CT-18 | Corrugated Tube- ID 22.0mm, L1.2m | |





Order Information

| Cat No. | Description |
|---------|------------------------------------|
| R419 | R419 Intelligent Animal Ventilator |

Included

| Cat No. | Description | Remarks |
|------------|--|----------|
| R419 | R419 Intelligent Animal Ventilator | |
| R409-TC | T-connector-9/16 | |
| R409-AH-80 | O ₂ hose-ID8.0-M16*1.5-9/16 connector, 1m | |
| RV-PEEP | Mechanical Ventilation PEEP Valve | Optional |

RWD Life Science | www.rwdstco.com

RWD SMD



Veterinary Multi-parameter Monitor

To monitor various parameters including SpO₂, ECG, NIBP, TEMP, RESP, ETCO₂, HR.





Embedded Motherboard

Long life & low maintenance cost





Digital Pulse Oximetry Accurate monitoring during low perfusion



Masimo Mainstream EtCO₂ (Optional) Reaction time <1S



Graphic Button and Backlight Design



Multiple Display Methods Meet individual needs



Built-in Recorder (Optional) Ready to print

Order Information

| Cat No. | Description | Remarks |
|---------|---|--------------------------------|
| RM300 | Veterinary Monitor with SpO₂, ECG, RESP, NIBP,TEMP and HR | |
| RM400M | Veterinary Monitor with EtCO ₂ , SpO ₂ , ECG、RESP、HR、NIBP and TEMP (Mainstream) | EtCO ₂ (Mainstream) |
| RM400S | Veterinary Monitor with EtCO ₂ , SpO ₂ , ECG、RESP、HR、NIBP and TEMP (Sidestream) | EtCO ₂ (Sidestream) |

R-CU500 Smart Animal ICU

R-CU500 The Smart Animal ICU, provides animal patients with the most comfortable medical care environment to help animals recover faster.



- Built-in preset parameters.
- Touchscreen LCD display for best simplicity.
- Air filter UVC lamp Anion function for effective sanitation management.
- Stable control of temperature, humidity, O2 concentration, and real-time CO2 concentration monitoring through electronic controller.
- Automatic ventilation based on CO2 concentration parameters and temperature; Built-in failure alarm covering abnormal high temperature, low humidity and CO2 levels.
- Support infusion / oxygen / atomization treatment, to meet various care conditions

Order Information

| Cat No. | Description |
|------------|------------------|
| R-CU500pro | Smart Animal ICU |

RM-C1 Veterinary Capnograph

To monitor ETCO2 and respiratory rate to help veterinarians identify animal ventilation in time.



- Real-time measurement of end-tidal carbon dioxide in animals
- Fast response, real-time EtCO2 waveform displayed in 5 seconds
- Audible and visual alarm system
- Clear display of waveform and value to provide doctors with comprehensive interpretation

Order Information

| Cat No. | Description |
|---------|-----------------------|
| RM-C1 | Veterinary Capnograph |

RWD Life Science | www.rwdstco.com www.rwdstco.com | RWD Life Science



ROC-5A Veterinary Oxygen Concentrator

ROC-5A oxygen concentrator is designed specifically for veterinary use. The high purity oxygen is separated from air by molecular sieve pressure swing adsorption(PSA)technology.



- Flow rate: 5L/min
- Generate medical oxygen (USP93%) from air, very cost effective
- With air filter, build-in filter and humidifier bottle, 99.999% of bacteria, viruses and odor can be removed effectively
- Low O₂ concentration alarm and abnormal alarm
- Quiet design to offer you a silent and comfortable working environment
- Suitable for supplying oxygen for anesthesia recovery, ICU, first aid etc

Order Information

| Cat No. | Description |
|---------|--------------------------------|
| ROC-5A | Veterinary Oxygen Concentrator |

Included

| Cat No. | Description | Remarks |
|--|-------------------------------------|---------|
| JAY-5BW | Oxygen Generator (Main Body) | |
| R510-25-86 | 9/16 TO Φ8 Cone Connector (1 piece) | |
| R-ST-60A Silicone Tubing- ID:6.0,OD:10.0mm | | |
| A5 | User Manual | |
| V-1 | Air Filter | |

Anesthesia Air Pump



Order Information

| Cat No. | Description |
|---------|--------------------------------|
| R510-29 | Anesthesia Air Pump (HF), 220V |
| R510-30 | Anesthesia Air Pump (HF), 110V |



Automatic Weighing Meter

Weighting RWD canister. Overweight alarm to avoid the waste gas pollution to your surgery room.



- Automatically weighing canister
- Alarm to remind of replacing canister
- One key to zero

Order Information

| Cat No. | Description |
|---------|---|
| R548 | Automatic Weighing Meter for Canister, American Standard* |
| R548-MB | Mounting Braket for Weighing Meter |

Gas filter Canister

Used to absorb isoflurane and sevoflurane. Net weight is 800g, and it can absorb 200g waste gas.



Order Information

| Cat No. | Description |
|------------|---|
| R510-31-6 | Gas Filter Canister, Large, Package of 6pcs |
| R510-31S-6 | Gas Filter Canister, Small, pkg of 6pcs |

13 RWD Life Science | www.rwdstco.com www.rwdstco.com | RWD Life Science | 14



▶ Vaporizer

RWD vaporizer adopts the most advanced design to accurately control the output concentration of anesthetic. Pour Fill, Easy Fill and Key Fill, Cagemont and Selectatec, Isoflurane and Sevoflurane available!







Pour-Fill

• Zero leakage, internal pressure up to 50kpa, safe and reliable

- >120ml capacity
- Isoflurane concentration range 0-5% (sevoflurane: 0-8%)
- Completely calibrating and testing one by one
- Built-in temperature compensator
- $_{\bullet}$ Working temperature range from 10 $^{\circ}\text{C}\,$ to 35 $^{\circ}\text{C}\,$





Order Information

| Cat No. | Description |
|---------|--|
| R580s | Vaporizer for Isoflurane - Pour Fill - Cagemount |
| R581S | Vaporizer for Sevoflurane- Pour Fill- Cagemount |
| R582S | Vaporizer for Sevoflurane - Easy Fill - Cagemount |
| R583S | Vaporizer for Isoflurane - Easy Fill - Cagemount |
| R580-IO | Vaporizer Endcaps (1 pair) |
| R580SS | Vaporizer for Isoflurane - Pour Fill - Selectatec |
| R581SS | Vaporizer for Sevoflurane - Pour Fill - Selectatec |
| R582SS | Vaporizer for Sevoflurane - Easy Fill - Selectatec |
| R583SS | Vaporizer for Isoflurane - Easy Fill - Selectatec |
| R580-SA | Vaporizer-suspension assembly-R640 |
| R584S | Vaporizer for Isoflurane- Key Fill-Cagemount |
| R585S | Vaporizer for Sevoflurane- Key Fill- Cagemount |
| R584SS | Vaporizer for Isoflurane - Key Fill - Selectatec |
| R585SS | Vaporizer for Sevoflurane- Key Fill- Selectatec |
| | |

▶ Anesthesia Induction Chambers







Order Information

| Cat No. | Description |
|---------|--------------------------------|
| V100-V | Anesthesia Induction Chamber-S |
| V101-V | Anesthesia Induction Chamber-M |
| V102-V | Anesthesia Induction Chamber-L |

▶ Shelf











R640-S1-MS-R

RAH-W

Order Information

| Cat No. | Description | Remarks |
|--------------|--|---------------|
| R-HS | Oxygen Concentrator Shelf (for H-base) | 340*484mm |
| R409-MB | Mounting Bracket for R409/R419 Ventilator | 235*220mm |
| R640-S1-MS-S | Monitor Shelf (for R640-S1) | 325*225mm |
| R640-S1-MS-R | Monitor Shelf (for R640-S1) | 325*225mm |
| RAH-W | Holder of veterinary anesthesia machine (Wall mount) | 265*150*120mm |

▶ Laryngoscopes

Apply to medium or large animals'laryngeal examination and tracheal intubation.

Order Information

| Cat No. | Description |
|-----------|--|
| C45001-00 | Miller Lamp Laryngoscope, 5 Blades+Handle+Case |
| C45002-00 | Miller F.O. Laryngoscope, 5 Blades+Handle+Case |





▶ Toggle Switch






R-GTS

R-TS R-GTS

▶ Regulator

Order Information





Toggle Switch for Oxygen Cylinder and Oxygen Concentrator

15 RWD Life Science | www.rwdstco.com | RWD Life Science | www.rwdstco.com | RWD Life Science | 16



▶ Mask


Apply to anesthesia induction, breathing rescue of felines, canines and other similar size animals.

Order Information

| Cat No. | Description |
|---------|--|
| 68637 | Small Feline Mask, 44mm(OD)*22mm(ID)*24mm(Depth) |
| 68638 | Medium Feline Mask, 57mm(OD)*30mm(ID)*29mm(Depth) |
| 68639 | Large Feline Mask, 60mm(OD)*33mm(ID)*35mm(Depth) |
| 68640 | Small Canine Mask, 87mm(OD)*30mm(ID)*73mm(Depth) |
| 68641 | Medium Canine Mask, 110mm(OD)*42mm(ID)*97mm(Depth) |
| 68642 | Large Canine Mask, 130mm(OD)*54mm(ID)*121mm(Depth) |
| 68643 | Feline/Canine Mask Set of 6 pcs |



▶ Breathing Bags



Order Information

| Cat No. | Description |
|--------------|--|
| R600-NL-0.5L | Reusable Non Latex Breathing Bag, 0.5L |
| R600-NL-1L | Reusable Non Latex Breathing Bag, 1L |
| R600-NL-2L | Reusable Non Latex Breathing Bag, 2L |
| R600-NL-3L | Reusable Non Latex Breathing Bag, 3L |

▶ Breathing Circuits



R-BC-DL











Order Information

| Cat No. | Description |
|-------------|---|
| R-BC-DL | Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m |
| R-NBC-F(JR) | Mapleson F Non-rebreathing Circuit (Jackson-Rees) |
| RSV-ET | Exhaust Tube |
| RSV-EK | Extension Kit |
| RSV-LDS | Low Dead Space Kit |
| RSV-GS | Gas Supply Hose |

► Air Hose







Order Information

| Cat No. | Description | Remarks | |
|-----------|---|------------------------------|-------------------|
| R-AH-80 | Air Hose (3-layer Compounded Rubber) - ID8.0-9/16(SGL)-2m | M16*1.5(M) and 9/16-18UNF(M) | C:- |
| R409-AH-1 | O ₂ Hose, Blue, 2*9/16 Connector, 1m | 9/16-18UNF(M)*2 | Green is optional |
| R409-AH-2 | O ₂ Hose, Blue, 2*9/16 Connector, 2m | 9/16-18UNF(M)*2 | |

► Avian Cole Endotubes



Order Information

| Cat No. | Description | Remarks |
|-------------|--------------------------------------|----------|
| R600-AET-10 | Silicone Avian Cole Endotube-ID1.0mm | 1pcs |
| R600-AET-15 | Silicone Avian Cole Endotube-ID1.5mm | 1pcs |
| R600-AET-20 | Silicone Avian Cole Endotube-ID2.0mm | 1pcs |
| R600-AET-25 | Silicone Avian Cole Endotube-ID2.5mm | 1pcs |
| R600-AET-30 | Silicone Avian Cole Endotube-ID3.0mm | 1pcs |
| R600-AET-35 | Silicone Avian Cole Endotube-ID3.5mm | 1pcs |
| R600-AET-40 | Silicone Avian Cole Endotube-ID4.0mm | 1pcs |
| R600-AET-50 | Silicone Avian Cole Endotube-ID5.0mm | 1pcs |
| R600-AET-60 | Silicone Avian Cole Endotube-ID6.0mm | 1pcs |
| R600-AET | Silicone Avian Cole Endotube(box) | 9pcs/box |

▶ Endotracheal Tubes



Order Information

| Cat No. | Description | Remarks |
|-------------|---|-----------|
| R600-ET | Cuffed Endotracheal Tubes, 17 sizes/Set | |
| R600-ET-20 | Endotracheal Tubes-I.D.2.0mm*O.D.3.3mm | 10pcs/box |
| R600-ET-25 | Cuffed Endotracheal Tubes-I.D.2.5mm*O.D.3.7mm | 10pcs/box |
| R600-ET-30 | Cuffed Endotracheal Tubes-I.D.3.0mm*O.D.4.0mm | 10pcs/box |
| R600-ET-35 | Cuffed Endotracheal Tubes-I.D.3.5mm*O.D.4.7mm | 10pcs/box |
| R600-ET-40 | Cuffed Endotracheal Tubes-I.D.4.0mm*O.D.5.3mm | 10pcs/box |
| R600-ET-45 | Cuffed Endotracheal Tubes-I.D.4.5mm*O.D.6.0mm | 10pcs/box |
| R600-ET-50 | Cuffed Endotracheal Tubes-I.D.5.0mm*O.D.6.7mm | 10pcs/box |
| R600-ET-55 | Cuffed Endotracheal Tubes-I.D.5.5mm*O.D.7.3mm | 10pcs/box |
| R600-ET-60 | Cuffed Endotracheal Tubes-I.D.6.0mm*O.D.8.0mm | 10pcs/box |
| R600-ET-65 | Cuffed Endotracheal Tubes-I.D.6.5mm*O.D.8.7mm | 10pcs/box |
| R600-ET-70 | Cuffed Endotracheal Tubes-I.D.7.0mm*O.D.9.3mm | 10pcs/box |
| R600-ET-75 | Cuffed Endotracheal Tubes-I.D.7.5mm*O.D.10.0mm | 10pcs/box |
| R600-ET-80 | Cuffed Endotracheal Tubes-I.D.8.0mm*O.D.10.7mm | 10pcs/box |
| R600-ET-85 | Cuffed Endotracheal Tubes-I.D.8.5mm*O.D.11.3mm | 10pcs/box |
| R600-ET-90 | Cuffed Endotracheal Tubes-I.D.9.0mm*O.D.12.0mm | 10pcs/box |
| R600-ET-95 | Cuffed Endotracheal Tubes-I.D.9.5mm*O.D.12.7mm | 10pcs/box |
| R600-ET-100 | Cuffed Endotracheal Tubes-I.D.10.0mm*O.D.13.3mm | 10pcs/box |
| R-ET-WHF | Wall Hanging Frame for Endotracheal Tubes (ID 2.0-13.0mm) | |
| | | |

▶ Urinary Catheter

/Cat Urinary Catheter/

Made of medical grade PVC material. Good biocompatibility, suitable for medium and long term indwelling (1-3 days).







Cat catheter without Stylet

Cat catheter with Stylet

Luer Lock

Order Information

| Product Name | Cat No. | Parameter(mm) | Remarks |
|--|--------------|---------------|---------|
| Cat Catheter Luer-1.0 (25pcs/box) | R600-CCD-10 | 1.0x130 | |
| Cat Catheter Luer-1.3 (25pcs/box) | R600-CCS-13 | 1.3x130 | |
| Cat Catheter-Open end-3Fr (10pcs/box) | R600-CCD1-10 | 1.0x130 | |
| Cat Catheter-Side hole-3Fr (10pcs/box) | R600-CCS1-10 | 1.0x130 | |

/Silicone Foley Catheter/

Made of 100% silicone material. High biocompatibility, suitable for long time indwelling (7 days). The catheter is visible under X-Ray.



Order Information

| Product Name | Cat No. | Parameter(mm) | Description |
|---|-------------|-----------------------------|-------------|
| Silicone Foley Catheter-2.0 (10pcs/box) | R600-FDC-20 | 6Fr (2.0x300)/ballon 3ml | |
| Silicone Foley Catheter-2.7 (10pcs/box) | R600-FDC-27 | 8Fr (2.7x300)/ballon 3-5ml | |
| Silicone Foley Catheter-3.3 (10pcs/box) | R600-FCD-33 | 10Fr (3.3x300)/ballon 3-5ml | |

/Dog Urinary Catheter/

Made of medical grade PVC material, good hardness, easy to intubate. The round tip and atraumatic side holes can protect the mucosal tissue during insertion.



Order Information

| Product Name | Cat No. | Parameter(mm) | Description |
|---|-------------|---------------|-------------|
| Disposable Dog Catheter-1.6 (10pcs/box) | R600-MDC-16 | 1.6x500 | |
| Disposable Dog Catheter-2.0 (10pcs/box) | R600-MDC-20 | 2.0x500 | |
| Disposable Dog Catheter-2.6 (10pcs/box) | R600-MDC-26 | 2.6x500 | |

Animal Surgical Instruments

SP0051-C Castration Surgery Instrument Kit for Pet

| 1 0051 C Castration Sargery mistrament Nictor 1 Ct | |
|---|---|
| Description | |
| Scalpel Handles with Ruler 3# Solid-12.5cm | 1 |
| 15# Scalpel Blades (Box of 100pcs) | 1 |
| 1x2 Teeth Tissue Forceps (Standard)-Str, Tip width 1.8mm, 14cm | 1 |
| Dressing Forceps-Str, Tip width 1.4mm, teeth length 19mm, 15.5cm | 1 |
| MAYO-STILLE Dissecting & Operating Scissors-B/B Str/37*9mm/14cm | 1 |
| MAYO-STILLE Dissecting & Operating Scissors-B/B Cvd/41*9.5mm/14cm | 1 |
| Operating Scissors w/Slender Blades-S/B Str/48.6*9.5mm/14.5cm | 1 |
| HALSTED Mosquito Forceps-Str, 21*1.2mm/12.5cm | 1 |
| HALSTED Mosquito Forceps-Cvd, 23.5*1mm/12.5cm | 1 |
| CRILE-RANKIN Artery Forceps-Str, 37*1mm/16cm | 1 |
| CRILE-RANKIN Artery Forceps-Cvd, 35*1.5mm/16cm | 1 |
| BACKHAUS Towel Forceps-11cm | 4 |
| Instrument Storage Portfolio, 32*22cm | 1 |
| Operating Scissors (Round Type)-S/S Str/48*10mm/15.5cm | 1 |
| MAYO-HEGAR Needle Holders-Str/21*2.8mm/16cm | 1 |
| ALLIS 4x5 Teeth Forceps-Str, 39*5mm Tips/15cm | 1 |



SP0050-S Sterilization Surgery Instrument Kit for Pets

| Description | |
|---|---|
| Scalpel Handles with Ruler 3# Solid-12.5cm | 1 |
| 15# Scalpel Blades (Box of 100pcs) | 1 |
| 1x2 Teeth Tissue Forceps (Standard)-Str, Tip width 1.8mm, 14cm | 1 |
| Dressing Forceps-Str, Tip width 1.4mm, teeth length 19mm, 15.5cm | 1 |
| MAYO-STILLE Dissecting & Operating Scissors-B/B Str/37*9mm/14cm | 2 |
| MAYO-STILLE Dissecting & Operating Scissors-B/B Cvd/41*9.5mm/14cm | 1 |
| Operating Scissors (Round Type)-S/S Str/48*10mm/15.5cm | 1 |
| Operating Scissors w/Slender Blades-S/B Str/48.6*9.5mm/14.5cm | 1 |
| HALSTED Mosquito Forceps-Str, 21*1.2mm/12.5cm | 2 |
| HALSTED Mosquito Forceps-Cvd, 23.5*1mm/12.5cm | 2 |
| CRILE-RANKIN Artery Forceps-Str, 37*1mm/16cm | 2 |
| CRILE-RANKIN Artery Forceps-Cvd, 35*1.5mm/16cm | 2 |
| BACKHAUS Towel Forceps-11cm | 4 |
| Ovarian Hooks, 1 Ball Tooth, 45°Cvd,10mm depth, 16cm | 1 |
| ALLIS 4x5 Teeth Forceps-Str, 39*5mm Tips/15cm | 1 |
| MAYO-HEGAR Needle Holders-Str/21*2.8mm/16cm | 1 |
| Instrument Storage Portfolio, 32*22cm | 1 |



DS1005-06 Winged Elevator Set for Dog (Standard Handle), 6pcs

| Description | , |
|--|---|
| Winged Elevators with Standard Handle, 1mm | 1 |
| Winged Elevators with Standard Handle, 2mm | 1 |
| Winged Elevators with Standard Handle, 3mm | 1 |
| Winged Elevators with Standard Handle, 4mm | 1 |
| Winged Elevators with Standard Handle, 5mm | 1 |
| Winged Elevators with Standard Handle, 6mm | 1 |
| Instrument Storage Portfolio, 32*22cm | 1 |



DS1001-15 Instrument Set for Dog w/Standard Handle Elevators, 15pcs

| Description | |
|--|---|
| Winged Elevators with Standard Handle, 1mm | 1 |
| Winged Elevators with Standard Handle, 2mm | 1 |
| Winged Elevators with Standard Handle, 3mm | 1 |
| Winged Elevators with Standard Handle, 4mm | 1 |
| Winged Elevators with Standard Handle, 5mm | 1 |
| Elevators with Standard Handle, Straight Root Tip Pick | 1 |
| Extraction Forceps for Dog, Curved,12.5cm | 1 |
| Tartar Removing Forceps with Curved Beak, 11.5cm | 1 |
| Double Ended Fine Periosteal Elevator | 1 |
| Subgingival Curette for Dog | 1 |
| Jaquette Supragingival Scaler for Dog | 1 |
| Sickle Scaler for Dog | 1 |
| Probe/Shepherd-hook Explorer with 1-10mm Marker | 1 |
| Minnesota Retractor, 14cm | 1 |
| Instrument Storage Portfolio, 32*22cm | 1 |
| Mouth Gag - 8cm Opening | 1 |



DS1010-10 Instrument Set for Cat w/Standard Handle Elevators, 10pcs

| Description | |
|---|---|
| Straight Tip Elevators with Standard Handle, 1mm | 1 |
| Straight Tip Elevators with Standard Handle, 2mm | 1 |
| Probe/Shepherd-hook Explorer with 1-10mm Marker | 1 |
| Double Ended Fine Periosteal Elevator | 1 |
| Resorptive Lesion Probe for Cat | 1 |
| Sickle Scaler for Cat | 1 |
| Curette for Small Dog or Cat | 1 |
| Mini Extraction Forceps for Cat, Straight, 12.5cm | 1 |
| Small Dog or Cat Right Angle Extraction Forceps, 12.5cm | 1 |
| Mouth Gag - 5cm Opening | 1 |
| Instrument Storage Portfolio, 32*22cm | 1 |



DS1006-06 Winged Elevator Set for Dog (Stubby Handle), 6pcs

| Description | |
|--|---|
| Winged Elevators with Stubby Handle, 1mm* | 1 |
| Winged Elevators with Stubby Handle, 2mm* | 1 |
| Winged Elevators with Stubby Handle, 3mm* | 1 |
| Winged Elevators with Stubby Handle, 4mm* | 1 |
| Winged Elevators with Stubby Handle, 5mm* | 1 |
| Winged Elevators with Stubby Handle, 6mm* | 1 |
| Instrument Storage Portfolio, 32*22cm | 1 |
| UK-SC002 Instrument sterilization box (Optional) | 1 |



RWD RWD


Veterinary Dental One-stop Solution

- Applicable to animal dental clinics
- One-stop solution for anesthesia system, dental unit and dental X-ray



Veterinary Anesthsia One-stop Solution

- Applicable to pet hospitals/clinics
 One-stop solution for anesthesia, monitoring, mechanical ventilation and waste gas scavenging
- Compact design and space saving
 Basket purchase, low maintenance cost



23 RWD Life Science | ww.rwdstco.com www.rwdstco.com | RWD Life Science