


Desktop Digital Stereotaxic Instruments

Digital displacement display, higher accuracy, higher efficiency!



Composition:

68025 and 68026 are complete sets of stereotaxic (see figure), both containing 68861N base plate, 68401 three-dimensional manipulator - left, 68021 rat adaptor, 68301 rat 18 ° ear bar, 68201 standard holder, 68026 includes an extra 68402 three-dimensional manipulator-right and 68023 digital display module.

Order Information:

Model	Product Description	Remarks
68025	Stereotaxic for Rat, SGL M, Digital, 18 Deg	
68026	Stereotaxic for Rat, Dual M, Digital, 18 Deg	
68027	Stereotaxic for Rat, SGL M, Digital, 45 Deg	
68028	Stereotaxic for Rat, Dual M, Digital, 45 Deg	
68018	Stereotaxic for Mouse, SGL M, Digital, 68030 Adaptor Incl	
68019	Stereotaxic for Mouse, Dual M, Digital, 68030 Adaptor Incl	
68045	Stereotaxic for Mouse, SGL M, Digital, 68055 Adaptor, 60 Deg	
68046	Stereotaxic for Mouse, Dual M, Digital, 68055 Adaptor, 60 Deg	
68537	Stereotaxic for Mouse, SGL M, Digital, 68077 Adaptor, 60 Deg	
68538	Stereotaxic for Rat, SGL M, Digital, 68078 Adaptor, 18 Deg	
68861N	Stereotaxic Base Plate, Standard	
68023	3-axis Digital Readout for Upgrade-68000	With digital base and displacement sensor