



## Halogen Fiber-lite Illuminator

The halogen fiber-lite illuminator adopts dual-piple pipe optic coaxial lighting, provide better field of view, with good brightness; stable output voltage, safe and reliable light bulb; size is about 280 x 100 x 165mm, and the light intensity can be adjusted continuously , with color temperature 3200K.

### Features

- ❑ Double optical fiber coaxial illumination, good field of view brightness.
- ❑ The output voltage is stable, and the bulb is safe and reliable to use.
- ❑ The appearance is small and beautiful, the size is about 280X100X165mm.
- ❑ Product use: animal surgery lighting, microscope lighting, etc.
- ❑ It can be continuously dimmed, and the color temperature is 3200K.
- ❑ Bulb: 21V 150W halogen lamp, rated life is about 200h.
- ❑ Cooling method: air cooling.
- ❑ Weight is about 3.0kg.

### Order Information

Model	Product Description	Remarks
76301	Fiber-Lite Illuminator(150W)	Dual-Pipe Fiber Optic , Bulb and Main Engine
76304	Main Engine	
76310	Dual-Pipe Fiber Optic ( Included in illuminator )	
76306	Replacement Bulb for 76310 Illuminator ( 21V 150W ) - OSRAM	