


R520 Portable Anesthesia Machine

Regular assembly (R540-A), vaporizer(R580), Inlet/Outlet End Caps(R580-IO), induction chamber, fixed bracket for gas evacuation apparatus, air pump, and portable assembly(R510-S), tubing are included in the R520 anesthesia machine.

Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatus should be purchased separately.



Order Information:

Cat No.	Product Description	Configuration	Remarks
R520IE	Movable Small Animal Anesthesia Machine-Isoflurane,Easy Fil	●	Including fixed bracket for induction chamber, evacuation apparatus, filter canister.
R520IP	Movable Small Animal Anesthesia Machine-Isoflurane,Pour Fil	●	
R520SE	Movable Small Animal Anesthesia Machine-Sevoflurane,Easy Fil	●	
R520SP	Movable Small Animal Anesthesia Machine-Sevoflurane,Pour Fil	●	
R600-N2O-4E	Expansion Assembly with 0.1-4LPM N ₂ O Flowmeter	○	
R510-09	2ch Flowmeter for Anesthesia, range 0.1-1.0L/min each	○	simultaneous anesthesia for two animals (adjust independently)

Note: ●Recommended configuration ○ Optional