


R550 Multi-output Anesthesia Machine

Regular assembly(R550-A), vaporizer(R580), Inlet/Outlet End Caps(R580-IO) and tube accessory kit are included in the R550 anesthesia machine. Different air supply, induction chambers, small animal masks, filter canisters and gas evacuation apparatuses can be selected to form the experimental platform for anesthesia recovery of different animals.

With the latest design, it can be applies to simultaneous anesthesia of multiple mice, rats, cats, rabbits and other animals (1-5 small animals), and every channel can be controlled independently



Order Information:

Cat No.	Product Description	Configuration	Remarks
R550IE	5ch Small Animal Anesthesia Machine-Isoflurane,Easy Fil	●	Contains connection tubing for mask & induction chamber
R550IP	5ch Small Animal Anesthesia Device-Isoflurane,Pour Fil	●	
R550SE	5ch Small Animal Anesthesia Machine-Sevoflurane,Easy Fil	●	
R550SP	5ch Small Animal Anesthesia Machine-Sevoflurane,Pour Fil	●	
R510-S	Portable Stand for R510 Series Anesthesia Machine	○	Upgrade to portable type
R600-N2O-4E	Expansion Assembly with 0.1-4LPM N ₂ O Flowmeter	○	

Note: ●Recommended configuration ○ Optional