

R620-MRI MR Conditional Animal Anesthesia Machine

R620-MRI has been tested by Shenzhen Institutes of Advanced Technology. And R620-MRI is proved not to produce displacement and not affect the operation and imaging in MRI systems of 3 Tesla and below 3 Tesla.



- Can be used in MRI systems of 3 Tesla and below 3 Tesla.
- High accuracy vaporizer.
- Good waste gas scavenging system.
- Complete function with both breathing and non-rebreathing circuits.

Order Information

Cat No.	Description
R620-MRI	MR Conditional Animal Anesthesia Machine

Included

Cat No.	Description	Remarks
R580S-MRI	Vaporizer for Isoflurane- Pour Fill- Cagemount-MRI	
R580-IO-MRI	Vaporizer Endcaps (1 pair)- MRI	
R510-31	Gas Filter Canister, Large	
R600-ET	Cuffed Endotracheal Tubes (1 set)	
R600-NL-0.5L	Reusable Non Latex Breathing Bag, 0.5L	
R600-NL-1L	Reusable Non Latex Breathing Bag, 1L	
R600-NL-2L	Reusable Non Latex Breathing Bag, 2L	
R-BC-DL	Breathing Circuit Duo-limb, 22mm(M)/15mm(F), 1.5m	
R-NBC-F(JR)	Mapleson F Non-rebreathing Circuit (Jackson-Rees)	
R-AH-MRI	Air Hose, ID8.0-9/16(SGL), 5m, MRI-Compatible	
R-CT-18	Corrugated Tube- ID 22.0mm, L1.2m	