

R640 Flexible Veterinary Anesthesia Machine-Wall Mount

R640 is an excellent choice for Non-operating room anesthesia, like Treatment Room, Imaging Room, Anesthesia Preparation Room and Feline Specialist Clinics. It can effectively improve the treatment rate and the medical experience of the owner and animals.



- Wall mount, ideal if you have limited space.
- Non-rebreathing circuit, ensure fast induction, quicker wake up and recovery.
- 10L flow meter, good for animals less than 50kg.
- Oxygen flush with safety restrictor.

Order Information

Cat No.	Description
R640	Flexible Veterinary Anesthesia Machine-Wall Mount

Included

Cat No.	Description	Remarks
R640-H	R640 Series Mounting Assembly	
R582S	Vaporizer for Sevoflurane- Easy Fill- Cagemount	Can be changed
R580-IO	Vaporizer Endcaps (1 pair)	
R-NBC-F(JR)	Mapleson F Non-rebreathing Circuit (Jackson-Rees)	
R600-NL-3L	Reusable Non Latex Breathing Bag, 3L	
R600-NL-1L	Reusable Non Latex Breathing Bag, 1L	
R-CT-18	Corrugated Tube- ID 22.0mm, L1.2m	