

RWD



 **Simplify Operation**
R640-S1
Veterinary Anesthesia Machine

RWD Vet Anesthesia, Easy & Safe

R640-S1 Veterinary Anesthesia Machine

Reliable Quality & Friendly Price

- **Compact structure & simple operation** — Simplify the veterinarian's working process.
- **Support Re-breathing circuit and Non-rebreathing circuit** (Basic configuration Re-breathing circuit 1.5m).
- **Pressure gauge** — Constantly monitoring the circuit pressure.
- **Recessed oxygen flush button** — Preventing emergencies.
- **Red rotary float** — More eye-catching and readable.
- **5 Freewheel with lock** — More stable.
- **3 years warranty** on main parts.



Configuration

Item	Specification
Vaporizer	Out concentration:0~5%(Isoflurane);0~8(Sevoflurane) Work temperature:10~35°C Volume>120ML
Flowmeter	0.1-4L/min
Breathing bag	1L
Re-breathing circuit	Low air resistance limbo breathing circuit,1.5m



Cylindrical Float
Stable Easy to read



CO₂ Absorbent
Usage Indicator



Embedded O₂
Flush Button



Rail-mounted CO₂
absorber

RWD Life Science Co.,Ltd

Add: 19-20/F, Building 7A, Shenzhen International Innovation Valley, Dashi 1st Road, Nanshan District, Shenzhen, Guangdong, P.R. China **E-mail:** rwd@rwdstco.com

RWD Life Science Inc.

Add: 850 New Burton Road, Suite 201, Dover, DE 19904 (New Office)
Tel: (858)900-5879 **Support:** service@rwdls.com