


Standard Stereotaxic Instruments

Classical U-type base design, professional services in the laboratory for 17 years!



Composition:

68001 and 68002 are complete sets of stereotaxic (see figure), both containing 68861N base plate, 68401 three-dimensional manipulator - left, 68021 rat adaptor, 68301 rat 18° ear bar, 68201 standard holder, and 68002 includes an extra 68402 three-dimensional manipulator - right.

Order Information:

Model	Product Description	Remarks
68001	Stereotaxic for Rat, SGL M, 18 Deg	
68002	Stereotaxic for Rat, Dual M, 18 Deg	
68005	Stereotaxic for Rat, SGL M, 45 Deg	
68006	Stereotaxic for Rat, Dual M, 45 Deg	
68037	Stereotaxic for Mouse, SGL M, 68030 Adaptor Incl	
68038	Stereotaxic for Mouse, Dual M, 68030 Adaptor Incl	
68043	Stereotaxic for Mouse, SGL M, 68055 Adaptor, 60 Deg	
68044	Stereotaxic for Mouse, Dual M, 68055 Adaptor, 60 Deg	
68535	Stereotaxic for Mouse, SGL M, 68077 Adaptor, 60 Deg	
68536	Stereotaxic for Rat, SGL M, 68078 Adaptor, 18 Deg	
68861N	Stereotaxic Base Plate, Standard	